

claricuff™

Blood Pressure Cuffs

Claricuff™ is an unrivaled level of exceptional product quality. The Neonatal Cuffs are especially made with safe materials and manufactured carefully for the health and well-being of Newborns and Babies.

Features

- Claricuff™ Neonatal Blood Pressure Cuffs are made for single patient use (disposable style)
- Claricuff™ Neonatal Blood Pressure Cuffs are made to AHA/AAMI Standards
- Claricuff™ Neonatal Blood Pressure Cuffs have different colors to aid in sizing
- Claricuff™ Neonatal Blood Pressure Cuffs are available in all standard sizes
- Claricuff™ Neonatal Blood Pressure Cuffs are made with highest quality materials
- Claricuff™ Neonatal Blood Pressure Cuffs do not contain natural latex rubber, DEHP or PVC



Claricuff™ Neonatal

| Description | Singe Tube | Double Tube | Cuff Size |
|---|------------|-------------|-----------|
| Claricuff™ BP Cuff TPU Neonatal #1 with ML fitting, 10/BG | CNN1ST-10 | CNN1DT-10 | size 1 |
| Claricuff™ BP Cuff TPU Neonatal #2 with ML fitting, 10/BG | CNN2ST-10 | CNN2DT-10 | size 2 |
| Claricuff™ BP Cuff TPU Neonatal #3 with ML fitting, 10/BG | CNN3ST-10 | CNN3DT-10 | size 3 |
| Claricuff™ BP Cuff TPU Neonatal #4 with ML fitting, 10/BG | CNN4ST-10 | CNN4DT-10 | size 4 |
| Claricuff™ BP Cuff TPU Neonatal #5 with ML fitting, 10/BG | CNN5ST-10 | CNN5DT-10 | size 5 |

Cuffs come standard with the ML fitting. It is offered with single tube or double tube.



Tapmedic LLC
4600 Touchton Road, Suite 1150
Jacksonville, FL 32246 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com