

# claricuff™

## Blood Pressure Cuffs

**Claricuff™ is an unrivaled level of exceptional product quality. The Neonatal Cuffs are especially made with safe materials and manufactured carefully for the health and well-being of Newborns and Babies.**

### Features

- Claricuff™ Neonatal Blood Pressure Cuffs are made for single patient use (disposable style)
- Claricuff™ Neonatal Blood Pressure Cuffs are made to AHA/AAMI Standards
- Claricuff™ Neonatal Blood Pressure Cuffs have different colors to aid in sizing
- Claricuff™ Neonatal Blood Pressure Cuffs are available in all standard sizes
- Claricuff™ Neonatal Blood Pressure Cuffs are made with highest quality materials
- Claricuff™ Neonatal Blood Pressure Cuffs do not contain natural latex rubber, DEHP or PVC



## Claricuff™ Neonatal

Description	Singe Tube	Double Tube	Cuff Size
Claricuff™ BP Cuff TPU Neonatal #1 with ML fitting, 10/BG	CNN1ST-10	CNN1DT-10	size 1
Claricuff™ BP Cuff TPU Neonatal #2 with ML fitting, 10/BG	CNN2ST-10	CNN2DT-10	size 2
Claricuff™ BP Cuff TPU Neonatal #3 with ML fitting, 10/BG	CNN3ST-10	CNN3DT-10	size 3
Claricuff™ BP Cuff TPU Neonatal #4 with ML fitting, 10/BG	CNN4ST-10	CNN4DT-10	size 4
Claricuff™ BP Cuff TPU Neonatal #5 with ML fitting, 10/BG	CNN5ST-10	CNN5DT-10	size 5

Cuffs come standard with the ML fitting. It is offered with single tube or double tube.



Tapmedic LLC  
9424 Baymeadows Rd, Suite 250  
Jacksonville, FL 32256 USA  
T +1 904-351-0244 F +1 904-351-0245  
www.tapmedic.com | info@tapmedic.com