



clariflo™
Pressure Infuser



Clariflo™ **Classic/ Ultimate/ Flex**

Designed for rapid IV infusion and pressure monitoring procedures.

Tapmedic is a global provider of high-quality medical equipment and healthcare devices, with a broad range of products used for rapid infusion, NIBP diagnostics, blood pressure cuffs and monitoring, and patient care.



Tapmedic LLC
9428 Baymeadows Rd, Suite 502
Jacksonville, FL 32256 USA
T +1 904-351-0244
www.tapmedic.com | info@tapmedic.com

ADVANTAGES

Quality you can count on

- No natural rubber latex
- CE Mark
- Durable, state of the art nylon bag
- Leak proof, premium seams for tear resistance
- Product needs minimal operator training
- 100% leak tested after manufacturing

Product features for easy use

- Integrated fluid bag hook
- Large hand bulb for quick inflation
- Security clamp
- IV pole loop with extended suspension
- No-crimp tubing

Choices: we have what you want

- 2 lines of products: Clariflo™ and Clariflo™ Flex
- 3 bag sizes: 500ml, 1000ml, 3000ml
- Disposable / reusable versions
- Piston gauge or aneroid gauge
- Stopcock or thumbwheel valve



Reinforced **see-through** TPU sleeve allows proper label observation.

OPTIONS

Pressure Gauges

The Clariflo™ is available in 2 versions of pressure gauges:



PISTON GAUGE

- Pressure range 0-300mmHg
- 360° gauge visibility
- Bleed valve prevents over-inflation



ANEROID GAUGE

- Pressure range 0-300mmHg
- Large, easy to read display
- Exact pressure control

Pressure Control Valves

The Clariflo™ offers a choice of 2 pressure control valves:



3-WAY STOPCOCK

- Precise control of pressure
- Conveniently operated by one hand



THUMBWHEEL VALVE


- Instant control of pressure
- Open or close mode

Clariflo™ Classic

The Clariflo™ Classic features a durable white nylon bag. It is made for single patient use and completely disposable, eliminating the need to clean or repair.

disposable	Part Number	Description	Box Qty
	500PGS	Clariflo™ Classic - 500ml bag with piston gauge and stopcock	2x10/box
	1000PGS	Clariflo™ Classic - 1000ml bag with piston gauge and stopcock	2x10/box
	3000PGS	Clariflo™ Classic - 3000ml bag with piston gauge and stopcock	2x5/box
	500PGT	Clariflo™ Classic - 500ml bag with piston gauge and thumbwheel valve	2x10/box
	1000PGT	Clariflo™ Classic - 1000ml bag with piston gauge and thumbwheel valve	2x10/box
	3000PGT	Clariflo™ Classic - 3000ml bag with piston gauge and thumbwheel valve	2x5/box
	500AGS	Clariflo™ Classic - 500ml bag with aneroid gauge and stopcock	2x10/box
	1000AGS	Clariflo™ Classic - 1000ml bag with aneroid gauge and stopcock	2x10/box
	3000AGS	Clariflo™ Classic - 3000ml bag with aneroid gauge and stopcock	2x5/box
	500AGT	Clariflo™ Classic - 500ml bag with aneroid gauge and thumbwheel valve	2x10/box
	1000AGT	Clariflo™ Classic - 1000ml bag with aneroid gauge and thumbwheel valve	2x10/box
	3000AGT	Clariflo™ Classic - 3000ml bag with aneroid gauge and thumbwheel valve	2x5/box

➤ 3 bag sizes available




Clariflo™ Ultimate

The Clariflo™ Ultimate is designed for reuse, easy cleaning and made with a durable, reinforced blue nylon bag.

reusable	Part Number	Description	Box Qty
	RU500PGS	Clariflo™ Ultimate - 500ml bag with piston gauge and stopcock	2x10/box
	RU1000PGS	Clariflo™ Ultimate - 1000ml bag with piston gauge and stopcock	2x10/box
	RU3000PGS	Clariflo™ Ultimate - 3000ml bag with piston gauge and stopcock	2x5/box
	RU500PGT	Clariflo™ Ultimate - 500ml bag with piston gauge and thumbwheel valve	2x10/box
	RU1000PGT	Clariflo™ Ultimate - 1000ml bag with piston gauge and thumbwheel valve	2x10/box
	RU3000PGT	Clariflo™ Ultimate - 3000ml bag with piston gauge and thumbwheel valve	2x5/box
	RU500AGS	Clariflo™ Ultimate - 500ml bag with aneroid gauge and stopcock	2x10/box
	RU1000AGS	Clariflo™ Ultimate - 1000ml bag with aneroid gauge and stopcock	2x10/box
	RU3000AGS	Clariflo™ Ultimate - 3000ml bag with aneroid gauge and stopcock	2x5/box
	RU500AGT	Clariflo™ Ultimate - 500ml bag with aneroid gauge and thumbwheel valve	2x10/box
	RU1000AGT	Clariflo™ Ultimate - 1000ml bag with aneroid gauge and thumbwheel valve	2x10/box
	RU3000AGT	Clariflo™ Ultimate - 3000ml bag with aneroid gauge and thumbwheel valve	2x5/box

➤ 3 bag sizes available



Clariflo™ Flex

The newest line is the Flex. Clariflo™ Flex provides the separate cuff and inflator. The concept is - use and keep the inflator. The cuff can be exchanged. It is semi-reusable and a reduction in bioburden.

4 styles available:

- Semi-disposable with piston gauge
- Semi-disposable with aneroid gauge
- Semi-reusable with piston gauge
- Semi-reusable with aneroid gauge

➤ **3 bag sizes available**



	Part Number	Description	Box Qty
disposable	500PGSFlex	Clariflo™ Flex Classic - 500ml bag with piston gauge and stopcock	2x10/box
	1000PGSFlex	Clariflo™ Flex Classic - 1000ml bag with piston gauge and stopcock	2x10/box
	3000PGSFlex	Clariflo™ Flex Classic - 3000ml bag with piston gauge and stopcock	2x5/box
	500AGSFlex	Clariflo™ Flex Classic - 500ml bag with aneroid gauge and stopcock	2x10/box
	1000AGSFlex	Clariflo™ Flex Classic - 1000ml bag with aneroid gauge and stopcock	2x10/box
	3000AGSFlex	Clariflo™ Flex Classic - 3000ml bag with aneroid gauge and stopcock	2x5/box
reusable	RU500PGSFlex	Clariflo™ Flex Ultimate - 500ml bag with piston gauge and stopcock	2x10/box
	RU1000PGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with piston gauge and stopcock	2x10/box
	RU3000PGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with piston gauge and stopcock	2x5/box
	RU500AGSFlex	Clariflo™ Flex Ultimate - 500ml bag with aneroid gauge and stopcock	2x10/box
	RU1000AGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with aneroid gauge and stopcock	2x10/box
	RU3000AGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with aneroid gauge and stopcock	2x5/box

Clariflo™ Flex Components

	Part Number	Description	Box Qty
disposable	500Flex	Clariflo™ Flex Classic - 500ml bag	2x10/box
	1000Flex	Clariflo™ Flex Classic - 1000ml bag	2x10/box
	3000Flex	Clariflo™ Flex Classic - 3000ml bag	2x5/box
	PGSFlex	Clariflo™ Flex Classic - piston gauge and stopcock	4x5/box
	AGSFlex	Clariflo™ Flex Classic - aneroid gauge and stopcock	4x5/box
	reusable	RU500Flex	Clariflo™ Flex Ultimate - 500ml bag
RU1000Flex		Clariflo™ Flex Ultimate - 1000ml bag	2x10/box
RU3000Flex		Clariflo™ Flex Ultimate - 3000ml bag	2x5/box
RUPGSFlex		Clariflo™ Flex Ultimate - piston gauge and stopcock	4x5/box
RUAGSFlex		Clariflo™ Flex Ultimate - aneroid gauge and stopcock	4x5/box

Tapmedic reserves the right to make changes to this sheet and the product at any time without notice. Clariflo™ is a registered trademark of Tapmedic.