



clariflo™ lite

Pressure Infusion Bag

Clariflo™ Lite

Designed for rapid IV infusion and pressure monitoring procedures.

Tapmedic is a global provider of high-quality medical equipment and healthcare devices, with a broad range of products used for rapid infusion, NIBP diagnostics, blood pressure cuffs and monitoring, and patient care.



Tapmedic LLC
9428 Baymeadows Rd, Suite 502
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com

clariflo™ lite

Pressure Infusion Bag

Clariflo™ Lite Features

ADVANTAGES

Quality you can count on

- No natural rubber latex
- CE Mark
- Durable, state of the art nylon bag
- Leak proof, premium seams for tear resistance
- Product needs minimal operator training
- 100% leak tested after manufacturing

Product features for easy use

- Integrated fluid bag hook
- Large hand bulb for quick inflation
- IV pole loop with extended suspension
- No-crimp tubing



PISTON GAUGE

- Pressure range 0-300mmHg
- 360° gauge visibility
- Bleed valve prevents over-inflation



3-WAY STOPCOCK

- Precise control of pressure
- Conveniently operated by one hand

Introducing the Clariflo™ Lite: Simple, Economical, Effective

The Clariflo™ Lite is designed to deliver the essentials, making it a more cost-effective option for hospitals and clinics. Featuring a durable white nylon bag, it offers reliable performance while remaining lightweight and easy to handle. Engineered for single-patient use, the Clariflo™ Lite is fully disposable, eliminating the need for cleaning or repairs—saving both time and resources.

With only the most essential features included, the Clariflo™ Lite provides an efficient, affordable solution without compromising on quality.

Part Number	Description	Box Qty
500PGSLT	Clariflo™ Lite - 500ml bag with piston gauge and stopcock	2x10/box
1000PGSLT	Clariflo™ Lite - 1000ml bag with piston gauge and stopcock	2x10/box