



ADAPTORS AND FITTINGS FOR THE BLOOD PRESSURE CUFFS CONNECTIONS CAN BE A REAL PUZZLE.



ADAPTORS AND FITTINGS


TAPMEDIC HAS SOLUTIONS AND CAN ASSIST IN THIS TASK.



Patient monitor hoses can be modified or adapters used to convert the patient monitoring equipment for use with the appropriate blood pressure cuff. Y adapters can be used to convert dual tube hoses to be used with a single tube BP cuff. Fittings can be used to for direct connections from monitors to the blood pressure cuffs.




www.tapmedic.com



Y-Adaptor with 1/8" Barb Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-182YRX4BLK (Black)
- ADP-182YRX4BLU (Blue)
- ADP-182YRX4WHT (White)



Y-Adaptor with 5/32" Barb Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-5322YRX4BLK (Black)
- ADP-5322YRX4BLU (Blue)
- ADP-5322YRX4WHT (White)



Y-Adaptor with Dinamap Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-DM2YRX4BLK (Black)
- ADP-DM2YRX4BLU (Blue)
- ADP-DM2YRX4WHT (White)



Y-Adaptor with Marquette Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-MQ2YRX4BLK (Black)
- ADP-MQ2YRX4BLU (Blue)
- ADP-MQ2YRX4WHT (White)



Y-Adaptor with Male Marquette Fittings and Metal HP Connector - 1/BG

Available in:

- ADP-MQML2YRX4BLK (Black)
- ADP-MQML2YRX4BLU (Blue)
- ADP-MQML2YRX4WHT (White)



Y-Adaptor with Male-Female Marquette Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-MQMF2YRX4BLK (Black)
- ADP-MQMF2YRX4BLU (Blue)
- ADP-MQMF2YRX4WHT (White)



Y-Adaptor with Dinaclick Fittings and Metal HP Connector - 1/BG

Available in:


- ADP-DC2YRX4BLK (Black)
- ADP-DC2YRX4BLU (Blue)
- ADP-DCYRX4WHT (White)



Y-Connector - 1/BG

Available in:

- ADP-YTUBECLK (Black)
- ADP-YTUBECLU (Blue)
- ADP-YTUBEWHT (White)



Twin Tube only - 1/BG

Available in:

- ADP-TWINTUBECLK (Black)



Twin Tube with Dinamap Fittings and Male Marquette Fittings - 1/BG

Available in:

- ADP-DM2TMQMLBLK (Black)



Twin Tube with Dinamap Fittings and Male/Female Marquette Fittings - 1/BG

Available in:

- ADP-DM2TMQMLBLK (Black)



Twin Tube with Dinamap Fittings and Female Marquette Fittings - 1/BG


Available in:

- ADP-DM2TMQMLBLK (Black)



Twin Tube with Male Dinamap Fittings and Male Marquette Fittings - 1/BG

Available in:
● ADP-MDM2TMQMLBLK (Black)



Twin Tube with Male Dinamap Fittings and Male/Female Marquette Fittings - 1/BG

Available in:
● ADP-MDM2TMQMFBLK (Male/Female Marquette)



Metal HP Connector - 1/BG

Available in:
ZFT-QC8220 Connector 1/8" / 3mm Barb
ZFT-QC8230 Connector 5/32" / 4mm Barb
ZFT-QC8240S Connector 3/16" / 5mm Barb




Metal Male HP Fitting 3/16" Barb - 1/BG

Available in:
ZFT-QC8030




Metal Male Dinamap Fitting 5/32" Barb - 1/BG

Available in:
ZFT-CS901




Dina Click Fitting 5/32" Barb - 20/BG

Available in:
ZFT-DC-20




Dinamap Fitting 5/32" Barb - 20/BG

Available in:
ZFT-DM532-20




Female Locking Luer 5/32" Barb - 20/BG

Available in:
ZFT-FL532-20



HP Fitting 5/32" Barb - 20/BG

Available in:
ZFT-HP532-20



Male Slip Luer 5/32" Barb - 20/BG

Available in:
ZFT-ML532-20



Male Locking Luer 5/32" Barb - 20/BG

Available in:
ZFT-MLL532-20




Female Marquette Fitting 5/32" Barb - 20/BG

Available in:
ZFT-MQ532-20




Male Marquette Fitting 5/32" Barb - 20/BG

Available in:
ZFT-MQML-20



5/32" to 5/32" Barb Connector - 20/BG

Available in:
ZFT-532M532-20



1/8" to 5/32" Barb Connector - 20/BG

Available in:
ZFT-18M532-20

Tapmedic has a full selection for all blood pressure cuff fittings and adaptors for all your patient monitoring requirements.

The use of the straight and Y adaptors have been evaluated as an effective method of converting and standardizing patient monitors with blood pressure cuffs.

Please contact Tapmedic for all for blood pressure hose, cuff and adaptors needs. We can customize to meet any patient monitor hose application.

Information needed:

- Patient Monitor Hose – single or double tube and type of adaptor
- Blood Pressure Cuff - single or double tube and type of fitting



Tapmedic LLC
9424 Baymeadows Rd, Suite 250
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com